



Inx6 Antibody

Product Code	CSB-PA895180XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9VR82
Immunogen	Synthesized peptide derived from internal of Drosophila melanogaster (Fruit fly) Inx6 protein (124-137aa)
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Target Names	Inx6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.