



laeA Antibody

Product Code	CSB-PA934639XA01EBP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6TLK5
Immunogen	Recombinant Emericella nidulans (Aspergillus nidulans) laeA protein
Raised In	Rabbit
Species Reactivity	Emericella nidulans (Aspergillus nidulans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Emericella nidulans (Aspergillus nidulans)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	laeA