



Acp26Aa Antibody

Product Code	CSB-PA339376XA01DLT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33735
Immunogen	Recombinant Drosophila mauritiana (Fruit fly) Acp26Aa protein
Raised In	Rabbit
Species Reactivity	Drosophila mauritiana (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila mauritiana (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Acp26Aa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.