



Nelf-E Antibody

Product Code	CSB-PA853331XA01DMP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95ZE9
Immunogen	Recombinant Drosophila virilis (Fruit fly) Nelf-E protein
Raised In	Rabbit
Species Reactivity	Drosophila virilis (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila virilis (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Nelf-E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.