



Nmdar1 Antibody

Product Code	CSB-PA465084XA01DLK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B3P2E5
Immunogen	Recombinant Drosophila erecta (Fruit fly) Nmdar1 protein
Raised In	Rabbit
Species Reactivity	Drosophila erecta (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila erecta (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Nmdar1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.