



Mat89Ba Antibody

Product Code	CSB-PA015924XA01DMJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4QX57
Immunogen	Recombinant Drosophila simulans (Fruit fly) Mat89Ba protein
Raised In	Rabbit
Species Reactivity	Drosophila simulans (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila simulans (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Mat89Ba
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.