




pho Antibody

Product Code	CSB-PA848344XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8ST83
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) pho protein (1-250aa)
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Research Area	Others
Target Names	pho

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: pho antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.