



PRL Antibody

Product Code	CSB-PA018724XA01FQO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29234
Immunogen	Recombinant Neovison vison (American mink) (Mustela vison) PRL protein
Raised In	Rabbit
Species Reactivity	Neovison vison (American mink) (Mustela vison)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Neovison vison (American mink) (Mustela vison)
Lead Time	Made-to-order (14-16 weeks)
Target Names	PRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.