



I(2)37Cg Antibody

Product Code	CSB-PA896102XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9VIZ0
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) I(2)37Cg protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	I(2)37Cg
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.