



POLR2A Antibody

Product Code	CSB-PA018327XA01DXU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11414
Immunogen	Recombinant Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus) POLR2A protein
Raised In	Rabbit
Species Reactivity	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cricetulus griseus (Chinese hamster) (Cricetulus barabensis griseus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	POLR2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.