



RpS15Aa Antibody

Product Code	CSB-PA744316XA01DMR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XIM8
Immunogen	Recombinant Drosophila yakuba (Fruit fly) RpS15Aa protein
Raised In	Rabbit
Species Reactivity	Drosophila yakuba (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila yakuba (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RpS15Aa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.