



Sgf11 Antibody

Product Code	CSB-PA453842XA01DMR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4PJ01
Immunogen	Recombinant Drosophila yakuba (Fruit fly) Sgf11 protein
Raised In	Rabbit
Species Reactivity	Drosophila yakuba (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila yakuba (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Sgf11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.