



SCRB16 Antibody

Product Code	CSB-PA752519XA01BZL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Q6R1
Immunogen	Recombinant Anopheles gambiae (African malaria mosquito) SCR16 protein
Raised In	Rabbit
Species Reactivity	Anopheles gambiae (African malaria mosquito)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Anopheles gambiae (African malaria mosquito)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SCR16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.