



Sod1 Antibody

Product Code	CSB-PA886977XA01DMO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9U4X2
Immunogen	Recombinant Drosophila teissieri (Fruit fly) Sod1 protein
Raised In	Rabbit
Species Reactivity	Drosophila teissieri (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila teissieri (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Sod1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.