



CG4953 Antibody

Product Code	CSB-PA842578XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95TN1
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) CG4953 protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CG4953
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.