



I(3)72Ab Antibody

Product Code	CSB-PA892901XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9VUV9
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) I(3)72Ab protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	I(3)72Ab
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.