




TTYH2 Antibody

Product Code	CSB-PA871565XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSA4
Immunogen	Recombinant Homo sapiens (Human) TTYH2 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TTYH2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: TTYH2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 22.88 kDa

Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.