



Vm32E Antibody

Product Code	CSB-PA396222XA01DIF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4UM16
Immunogen	Recombinant Drosophila santomea (Fruit fly) Vm32E protein
Raised In	Rabbit
Species Reactivity	Drosophila santomea (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila santomea (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	Vm32E