



XDH Antibody

Product Code	CSB-PA026189XA01CAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08793
Immunogen	Recombinant Calliphora vicina (Blue blowfly) (Calliphora erythrocephala) XDH protein
Raised In	Rabbit
Species Reactivity	Calliphora vicina (Blue blowfly) (Calliphora erythrocephala)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Calliphora vicina (Blue blowfly) (Calliphora erythrocephala)
Lead Time	Made-to-order (14-16 weeks)
Target Names	XDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.