



Vm32E Antibody

Product Code	CSB-PA394559XA01DLT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4UM12
Immunogen	Recombinant Drosophila mauritiana (Fruit fly) Vm32E protein
Raised In	Rabbit
Species Reactivity	Drosophila mauritiana (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila mauritiana (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Vm32E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.