



GF15731 Antibody

Product Code	CSB-PA026480XA01DKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B3MPC0
Immunogen	Recombinant Drosophila ananassae (Fruit fly) GF15731 protein
Raised In	Rabbit
Species Reactivity	Drosophila ananassae (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila ananassae (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GF15731
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.