



aerA Antibody

Product Code	CSB-PA362609XA01AXW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09167
Immunogen	Recombinant Aeromonas hydrophila aerA protein
Raised In	Rabbit
Species Reactivity	Aeromonas hydrophila
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aeromonas hydrophila
Lead Time	Made-to-order (14-16 weeks)
Target Names	aerA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.