



txyA Antibody

Product Code	CSB-PA862754XA01VFF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9LCB9
Immunogen	Recombinant Vibrio sp. txyA protein
Raised In	Rabbit
Species Reactivity	Vibrio sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vibrio sp.
Lead Time	Made-to-order (14-16 weeks)
Target Names	txyA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.