



acpP Antibody

Product Code	CSB-PA015636XA01FQU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O34163
Immunogen	Recombinant Brachyspira hyodysenteriae (Treponema hyodysenteriae) acpP protein
Raised In	Rabbit
Species Reactivity	Brachyspira hyodysenteriae (Treponema hyodysenteriae)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Brachyspira hyodysenteriae (Treponema hyodysenteriae)
Lead Time	Made-to-order (14-16 weeks)
Target Names	acpP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.