



acpP Antibody

Product Code	CSB-PA015636XA01LFJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80920
Immunogen	Recombinant Leucothrix mucor acpP protein
Raised In	Rabbit
Species Reactivity	Leucothrix mucor
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Leucothrix mucor
Lead Time	Made-to-order (14-16 weeks)
Gene Names	acpP