



Protein 5NUC Antibody

Product Code	CSB-PA345587XA01FQY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52307
Immunogen	Recombinant Rhipicephalus microplus (Cattle tick) (Boophilus microplus) Protein 5NUC protein
Raised In	Rabbit
Species Reactivity	Rhipicephalus microplus (Cattle tick) (Boophilus microplus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhipicephalus microplus (Cattle tick) (Boophilus microplus)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.