



aacC3 Antibody

Product Code	CSB-PA680367XA01EZW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29808
Immunogen	Recombinant Pseudomonas aeruginosa aacC3 protein
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas aeruginosa
Lead Time	Made-to-order (14-16 weeks)
Target Names	aacC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.