



penDE Antibody

Product Code	CSB-PA325780XA01ETG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15802
Immunogen	Recombinant Penicillium chrysogenum (Penicillium notatum) penDE protein
Raised In	Rabbit
Species Reactivity	Penicillium chrysogenum (Penicillium notatum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Penicillium chrysogenum (Penicillium notatum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	penDE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.