



Aliphatic amidase Antibody

Product Code	CSB-PA874250XA01GFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RQ17
Immunogen	Recombinant Geobacillus stearothermophilus (Bacillus stearothermophilus) Aliphatic amidase protein
Raised In	Rabbit
Species Reactivity	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.