



aepA Antibody

Product Code	CSB-PA312376XA01ESP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06555
Immunogen	Recombinant Pectobacterium carotovorum subsp. carotovorum (Erwinia carotovora subsp. carotovora) aepA protein
Raised In	Rabbit
Species Reactivity	Pectobacterium carotovorum subsp. carotovorum (Erwinia carotovora subsp. carotovora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pectobacterium carotovorum subsp. carotovorum (Erwinia carotovora subsp. carotovora)
Lead Time	Made-to-order (14-16 weeks)
Target Names	aepA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.