



adh Antibody

Product Code	CSB-PA955927XA01SRL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50381
Immunogen	Recombinant Sulfolobus sp. (strain RC3) adh protein
Raised In	Rabbit
Species Reactivity	Sulfolobus sp. (strain RC3)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sulfolobus sp. (strain RC3)
Lead Time	Made-to-order (14-16 weeks)
Target Names	adh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.