



pcbAB Antibody

Product Code	CSB-PA311543XA01FOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01757
Immunogen	Recombinant Streptomyces clavuligerus pcbAB protein
Raised In	Rabbit
Species Reactivity	Streptomyces clavuligerus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces clavuligerus
Lead Time	Made-to-order (14-16 weeks)
Target Names	pcbAB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.