



BMY1 Antibody

Product Code	CSB-PA318056XA01INB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10537
Immunogen	Recombinant Ipomoea batatas (Sweet potato) (Convolvulus batatas) BMY1 protein
Raised In	Rabbit
Species Reactivity	Ipomoea batatas (Sweet potato) (Convolvulus batatas)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ipomoea batatas (Sweet potato) (Convolvulus batatas)
Lead Time	Made-to-order (14-16 weeks)
Target Names	BMY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.