



amyC Antibody

Product Code	CSB-PA334310XA01FRK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37729
Immunogen	Recombinant Thermoanaerobacterium thermosulfurigenes (Clostridium thermosulfurogenes) amyC protein
Raised In	Rabbit
Species Reactivity	Thermoanaerobacterium thermosulfurigenes (Clostridium thermosulfurogenes)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermoanaerobacterium thermosulfurigenes (Clostridium thermosulfurogenes)
Lead Time	Made-to-order (14-16 weeks)
Target Names	amyC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.