



cnbCa Antibody

Product Code	CSB-PA968839XA01CNN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6J1Z5
Immunogen	Recombinant Comamonas testosteroni (Pseudomonas testosteroni) cnbCa protein
Raised In	Rabbit
Species Reactivity	Comamonas testosteroni (Pseudomonas testosteroni)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Comamonas testosteroni (Pseudomonas testosteroni)
Lead Time	Made-to-order (14-16 weeks)
Target Names	cnbCa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.