



AMO Antibody

Product Code	CSB-PA319452XA01EVQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12807
Immunogen	Recombinant Pichia angusta (Yeast) (Hansenula polymorpha) AMO protein
Raised In	Rabbit
Species Reactivity	Pichia angusta (Yeast) (Hansenula polymorpha)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pichia angusta (Yeast) (Hansenula polymorpha)
Lead Time	Made-to-order (14-16 weeks)
Target Names	AMO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.