



# aglB Antibody

Product Code	CSB-PA530788XA01TNK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O86959
Immunogen	Recombinant Thermotoga neapolitana aglB protein
Raised In	Rabbit
Species Reactivity	Thermotoga neapolitana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermotoga neapolitana
Lead Time	Made-to-order (14-16 weeks)
Gene Names	aglB