



ahpA Antibody

Product Code	CSB-PA881692XA01ESG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9L8J3
Immunogen	Recombinant Pasteurella multocida (strain Pm70) ahpA protein
Raised In	Rabbit
Species Reactivity	Pasteurella multocida (strain Pm70)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pasteurella multocida (strain Pm70)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ahpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.