



amiE Antibody

Product Code	CSB-PA888675XA01BRE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9L543
Immunogen	Recombinant Bacillus sp. amiE protein
Raised In	Rabbit
Species Reactivity	Bacillus sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus sp.
Lead Time	Made-to-order (14-16 weeks)
Gene Names	amiE