



amiE Antibody

Product Code	CSB-PA164545XA01RLJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01360
Immunogen	Recombinant Rhodococcus erythropolis (Arthrobacter picolinophilus) amiE protein
Raised In	Rabbit
Species Reactivity	Rhodococcus erythropolis (Arthrobacter picolinophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodococcus erythropolis (Arthrobacter picolinophilus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	amiE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.