



apu Antibody

Product Code	CSB-PA334197XA01FRM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36905
Immunogen	Recombinant Thermoanaerobacterium saccharolyticum apu protein
Raised In	Rabbit
Species Reactivity	Thermoanaerobacterium saccharolyticum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermoanaerobacterium saccharolyticum
Lead Time	Made-to-order (14-16 weeks)
Target Names	apu
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.