



Aquaporin Antibody

Product Code	CSB-PA638811XA01FRN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q25074
Immunogen	Recombinant Haematobia irritans exigua (Buffalo fly) Aquaporin protein
Raised In	Rabbit
Species Reactivity	Haematobia irritans exigua (Buffalo fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Haematobia irritans exigua (Buffalo fly)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.