



arcC1 Antibody

Product Code	CSB-PA358297XA01ELX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A2X8
Immunogen	Recombinant Enterococcus faecium (Streptococcus faecium) arcC1 protein
Raised In	Rabbit
Species Reactivity	Enterococcus faecium (Streptococcus faecium)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterococcus faecium (Streptococcus faecium)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	arcC1