



aphA Antibody

Product Code	CSB-PA517483XA01SWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O08430
Immunogen	Recombinant Salmonella typhi aphA protein
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella typhi
Lead Time	Made-to-order (14-16 weeks)
Target Names	aphA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.