



NAT1 Antibody

Product Code	CSB-PA015464XA01MRG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50292
Immunogen	Recombinant Mesocricetus auratus (Golden hamster) NAT1 protein
Raised In	Rabbit
Species Reactivity	Mesocricetus auratus (Golden hamster)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mesocricetus auratus (Golden hamster)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	NAT1