



asr Antibody

Product Code	CSB-PA842390XA01EKS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93MH6
Immunogen	Recombinant Enterobacter cloacae asr protein
Raised In	Rabbit
Species Reactivity	Enterobacter cloacae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacter cloacae
Lead Time	Made-to-order (14-16 weeks)
Target Names	asr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.