



Atp5f1b Antibody

Product Code	CSB-PA002350XA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10719
Immunogen	Recombinant Rattus norvegicus (Rat) Atp5f1b protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus (Rat)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	Atp5f1b