




Bifunctional beta-D-glucosidase/beta-D-fucosidase Antibody

Product Code	CSB-PA308646XA01BGO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P94248
Immunogen	Recombinant Bifidobacterium breve Bifunctional beta-D-glucosidase/beta-D-fucosidase protein (253-460aa)
Species Reactivity	Bifidobacterium breve
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bifidobacterium breve
Target Names	N/A

Image



Western Blot

Positive WB detected in: recombinant protein
 All lanes: Bifunctional beta-D-glucosidase/beta-D-fucosidase Antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 47 kDa
 Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.