



caxA Antibody

Product Code	CSB-PA265480XA01FRT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0ZAI3
Immunogen	Recombinant Alkalimonas amylolytica caxA protein
Raised In	Rabbit
Species Reactivity	Alkalimonas amylolytica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Alkalimonas amylolytica
Lead Time	Made-to-order (14-16 weeks)
Target Names	caxA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.