



tpn39b Antibody

Product Code	CSB-PA329934XA01TPR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29725
Immunogen	Recombinant Treponema pallidum (strain Nichols) tpn39b protein
Raised In	Rabbit
Species Reactivity	Treponema pallidum (strain Nichols)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Treponema pallidum (strain Nichols)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tpn39b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.